



Oxygen and More...

INSPIRON™ PORTABLE OXYGEN CONCENTRATORS



INSPIRON™ – THE OBVIOUS CHOICE

- ▶ Smooth and continuous supply of high quality oxygen.
- ▶ Available in 5LPM and 8LPM variants.
- ▶ Optimal power consumption with standard electrical requirements.
- ▶ Quiet operation with highly controlled sound vibrations.
- ▶ Innovatively designed compact cabinet, made of aesthetic ABS plastic, mounted on castor wheels for easy movement.
- ▶ Cabinet dimensions make it possible to store the machine even in small areas.
- ▶ Rugged body of the product makes it easy for transportation.
- ▶ Controls are simple to read and operate.
- ▶ Panel indicator provides necessary information with clarity and effectiveness - Immediate display to indicate if oxygen output is as per the desired levels.



◀ INSPIRON™ 5 - 5L Oxygen Concentrator

Technical Specification

Product Weight	: 15.5 kg
Product Width	: 350mm
Product Length/Depth	: 340mm
Product Height	: 530mm
Oxygen Purity	: 93%±2%@5L,25°C
Flow Range(lpm)	: 0.5~5.0 lpm
Operating System	: Time Cycle/Pressure Swing
Sound Level	: <40dBA(ISO 8359:1966 From Front)
Electrical Requirements	: 220~240Volts AC,50Hz, 1.8Amp
Power Consumption	: 315 Watts Average

Quality and safety features

- ▶ Filters for dust , oil and moisture
- ▶ Thermal switch for high temperature cutoff
- ▶ Low purity AV alarm (optional)
- ▶ Oxygen purity indicator built in (optional)

Certifications

- ▶ ISO 9001 : 2008 Company
- ▶ ISO 13485 Manufacturing

Accessories Included

- ▶ Disposable Humidifier Bottle, Nasal Cannula, Operator Manual



▼ *inspiro*[™] 8 - 8L Oxygen Concentrator



Technical Specification

Product Weight	: 15.85 Kg
Product Width	: 355 mm
Product Length/Depth	: 305 mm
Product Height	: 508 mm
Oxygen Purity	: 93%@8L,25°C
Flow Range(lpm)	: 0.5~5.0 lpm
Operating System	: Time Cycle/Pressure Swing
Sound Level	: <48dBA(ISO 8359:1966 From Front)
Electrical Requirements	: 220Volts AC,50Hz,
Power Consumption	: 315 Watts Average

Optionals

- ▶ Oxygen Purity Indicator
- ▶ Flow Splitter
- ▶ Reusable Humidifier



helix
CORPORATION

Buy your *inspiro*[™] from:

878, 17th Cross, 9th 'B' Main Road
ISRO Layout, Bangalore - 560 078
Ph: 080-2686 0687 / 2686 0689
Fax: 080-2686 0692
E-mail: sales@inspiro2.com
http: inspiro2.com

I026112009 REV00